
Lime-Lite™ ENHANCED

BIOACTIVE LIGHT CURE CAVITY LINER

POWERED BY
CRYSTA™
MCP TECHNOLOGY
U.S. PATENT NO 11,116,701 B2

XP-LLE-IN-07w
REV: 10/2024

EN

PULPDENT®
LIME-LITE™ ENHANCED
Light Cure Cavity Liner
Calcium, Phosphate, and Fluoride Releasing
• Radiopaque • Moisture Tolerant • Universal
Dentin Shade • Contains no BPA
Lime-Lite Enhanced is a moisture tolerant, light cure, radiopaque cavity liner and base with self-adhesive properties. It is a mineral-enriched urethane dimethacrylate resin that releases calcium, phosphate, and fluoride ions and is formulated for use with adhesives, composites, and conventional restorative materials. It is similar in chemistry to dental composites and chemically bonds to composites and adhesive primers. When increased bond strength to dentin is desired, use of a bonding agent is recommended. Lime-Lite Enhanced contains a patented rubberized-resin component that absorbs shock and stress and resists fracture and chipping. The material sets hard upon light curing, has a dentin shade, and is virtually insoluble in water and oral fluids. Lime-Lite Enhanced contains no BHT, no Bisphenol A, no Bis-GMA, and no BPA derivatives.

PRODUCT DESCRIPTION
Lime-Lite Enhanced is a mineral-enriched urethane dimethacrylate resin cavity liner and base that is moisture friendly, light cure, and radiopaque. It releases calcium, phosphate, and fluoride.

INDICATIONS
Indicated as a light cure cavity liner and base.
CONTRAINDICATIONS
Not indicated for direct placement on the exposed pulp.

INSTRUCTIONS FOR USE
1. Prepare cavity, rinse, and dry, but do not desiccate the tooth.
2. Place a pink 20-gauge dispenser tip securely onto the Lime-Lite Enhanced syringe.
3. Apply 1-1½mm of Lime-Lite Enhanced to the dentin.
4. Light cure for 20 seconds (450 mW/cm2). Do not exceed manufacturer's recommended curing time for the light you are using.
5. If used as a base, apply additional increments of 1-1½mm and light cure each layer. Cap syringe immediately after use.

6. Proceed with restoration. Lime-Lite Enhanced can be etched and is compatible with all bonding agents and dental restorative materials. It can be used with both open sandwich and closed sandwich techniques.

NOTE:
• If using Lime-Lite Enhanced under water based cements, such as zinc phosphate, polyacrylic or glass ionomer cement, remove any uncured resin from the surface of Lime-Lite Enhanced (oxygen inhibited layer) with a pledget of cotton moistened with alcohol.
STORAGE AND HANDLING
• Store tightly sealed in original container at cool room temperature. Avoid direct light, extremes of temperature, contamination and sources of ignition.
• Shelf life of unopened product: 2 years from date of manufacture.
• Re-cap immediately after use.

CAUTION
Uncured material may cause eye or skin irritation on contact. Dental professionals should wear safety glasses and surgical gloves. Do not exceed manufacturer's recommended curing time for the light you are using.

REPORTING OF SERIOUS INCIDENTS
Any serious incidents involving the device should be reported to the manufacturer and the Competent Authority of the Member State in which the incident occurred, or other appropriate regulatory body. Reports of serious incidents should be sent to:
Pulpdent Corporation
80 Oakland Street, Watertown, MA 02472 USA
Tel: 617-926-6666 / Fax: 617-926-6262
email: pulpdent@pulpdent.com
Note: Apply disposable barrier sleeves/wraps over multiple-use dental dispensers before use with each patient. For additional information, refer to: <https://www.fda.gov/medical-devices/dental-devices/multiple-use-dental-dispenser-devices>
Note: Refer to the Safety Data Sheet for a list of materials present in this product.

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Ravelli S.p.A.
Hofmeester Dental B.V.
Lifco Dental AB
UNIDENT AB
PULPDENT® Corporation 80 Oakland Street • Watertown, MA 02472 • U.S.A.